

# Aspergillus GM

Lateral Flow Assay



REF: AF2003

CE: Serum, BAL

## It's About Time



**45 MIN  
to Results!**

### Test Features

- Results in Less than an Hour
- Reduces Time to Proper Patient Care
- Highly Sensitive & Specific




The LFA could aid in quick clinical decisions in ICU patients, which could lead to increased survival of patients with true IPA.<sup>5</sup>

# The IMMY sōna *Aspergillus* galactomannan LFA has been suggested as a diagnostic tool in the 2020 ECMM CAPA-IAA-ICU management algorithm guidelines.<sup>1</sup>

Good agreement to existing product, better agreement to EORTC (true disease status):

Performance of IMMY ASP LFA	Specimen	Sensitivity	Specificity	NPV
vs. Proven IA <sup>2</sup>	Serum & BAL	100%	95.5%	
vs. Probable IA <sup>2</sup>	Serum & BAL	87.5%	96.2%	
vs. Proven IPA <sup>3</sup>	BAL	91%	92%	99%
vs. EORTC/MSG <sup>4</sup>	Serum	91%	91%	
vs. 2008 EORTC/MSG <sup>5</sup>	BAL	88%	81%	94%
vs. EORTC excluding GM <sup>5</sup>	BAL	100%	81%	96%
vs. AspICU <sup>5</sup>	BAL	94%	81%	97%
vs. modified AspICU <sup>5</sup>	BAL	87%	81%	94%
vs. modified AspICU excluding GM <sup>5</sup>	BAL	97%	81%	98%

1. ECMM. COVID-19/Influenza-Associated Pulmonary Aspergillosis – Management. 2020  
 2. Jani K, McMillen T, Morjaria S, Babady NE. Performance of the sōna *Aspergillus* Galactomannan Lateral Flow Assay in a Cancer Patient Population. J Clin Microbiol. 2021; JCM0059821.  
 3. Mercier T, et al. Lateral flow assays for diagnosing invasive pulmonary aspergillosis in adult hematology patients: A comparative multicenter study. Med Mycol. 2019.  
 4. Serin I, Dogu MH. Serum *Aspergillus* galactomannan lateral flow assay for the diagnosis of invasive aspergillosis: A single-centre study. Mycoses. 2021;64(6):678-683.  
 5. Mercier T, et al. Point of care aspergillus testing in intensive care patients. Crit Care. 2020.



**LFA Cube Reader** provides objective numeric value and allows for easy batch testing.

## CONTACT US

Ready to evaluate today?  
 Contact us for a validation kit now.

Distributed by