

## Sorenson LTS Compatible Pipet Tips

Effortless ejection, superior sealing, ergonomic excellence.




The new Sorenson LTS Compatible Tips are ergonomically designed to reduce discomfort and optimized to work with Rainin LiteTouch™ pipettors, all with the quality you have come to expect from Sorenson.

■ **Comfort you can feel** : LTS Compatible Tips reduce the force required to seal and eject, enabling a streamlined experience from aspiration to dispensing. Avoid slamming with excessive force on the pipettor for a proper seal, eliminate the struggle to eject tips, and minimize hand and wrist strain after extensive pipetting

■ **Quality** : Our United States manufacturing facility is audited by third parties to ensure adherence to regulated manufacturing standards, including ISO 13485 and ISO 14001. Sterile products feature a Sterility Assurance Level (SAL) of  $10^{-3}$ , and all products are made using USP Class VI resin

■ **Convenience of choice** : LTS Compatible Tips are available in your favorite packaging formats, including racked, racked/filtered, tip station, reload stacks, and our new more sustainable solution: ClearLine® HybridRack™ Pipet Tips. If you're not using an LTS pipettor, no problem. We provide a complete range of universal pipet tips to choose from

■ **Peace of mind** : Our LTS Compatible pipet tips are RNase-/DNase-free, Human gDNA-free, PCR inhibition-free, and nonpyrogenic so you can pipet stress-free, knowing your samples are safe. Additionally, the tips are autoclavable at 121°C for 15 minutes



■ Inert polyethylene filter reduces the risk of crossover contamination, protecting your samples

■ Ultra Clear USP Class VI resin

■ Fine point delivery system ensures crisp sample dispensing

■ Flexible sealing area optimized to work with LTS pipettors enables less force to be used during pipetting

■ Graduation lines for quick accuracy checks

Distributed by

**BIOMEDICA**  
www.bmgrp.eu

## Ordering Information

### MultiRack, Filtered, Sterile



MultiRack

Cat. No.	Description	Qty/Pk	Qty/Cs
0087940	20 µL LTS compatible tip, filtered, racked, sterile	960	4,800
0087950	200 µL LTS compatible tip, filtered, racked, sterile	960	4,800
0087960	1000 µL LTS compatible tip, filtered, racked, sterile	768	3,840

### MultiRack, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
0087910	20 µL LTS compatible tip, racked, sterile	960	4,800
0087920	200 µL LTS compatible tip, racked, sterile	960	4,800
0087930	1000 µL LTS compatible tip, racked, sterile	768	3,840

### MultiRack

Cat. No.	Description	Qty/Pk	Qty/Cs
0087880	20 µL LTS compatible tip, racked	960	4,800
0087890	200 µL LTS compatible tip, racked	960	4,800
0087900	1000 µL LTS compatible tip, racked	768	3,840



HybridRack™

### Sorenson Reload Stack

Cat. No.	Description	Qty/Pk	Qty/Cs
0088000	20 µL LTS compatible tip, Sorenson Reload Stack	960	3840
0088010	200 µL LTS compatible tip, Sorenson Reload Stack	960	3840
0088020	1000 µL LTS compatible tip, Sorenson Reload Stack	768	3072

### Tip Station, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
0088040	20 µL LTS compatible tip, tip station, sterile	960	4,800
0088060	200 µL LTS compatible tip, tip station, sterile	960	4,800

### Tip Station

Cat. No.	Description	Qty/Pk	Qty/Cs
0088030	20 µL LTS compatible tip, tip station	960	4,800
0088050	200 µL LTS compatible tip, tip station	960	4,800

### HybridRack, Filtered, Sterile

Cat. No.	Description	Qty/Pk	Qty/Cs
0087970	20 µL LTS compatible tip, HybridRack, Filtered, Sterile	960	4800
0087980	200 µL LTS compatible tip, HybridRack, Filtered, Sterile	960	4800
0087990	1000 µL LTS compatible tip, HybridRack, Filtered, Sterile	768	3840



Reload Stack

Distributed by