

MicroCentrifuge Tubes "SafeSeal"



0,65 ml; 1,7 ml; 2 ml;

Superior clarity - glass like: easy viewing of sample and pellets

Certified as DNase-, RNase-, PCR inhibitor-, human DNA- and Pyrogen- free

High speed centrifugation: RCF rating of 24.000 x g

Temperature rating from -80 °C to 121 °C

Fit all standard rotors

Autoclavable

Positive seal design: **will not pop open during boiling or freezing**

External graduation marks

Frosted writing surfaces

Available in 6 colors including natural

Available low binding

2 ml tubes available in **Dolphin shape**:

ideal for pelletizing samples

RCF rating of 15.000 x g

