



Single- and Multi-channel Pipettors

- Volume locking prevents accidental counter shift
- Ergonomic : contoured handle to adapt to the hand, finger support to reduce wrist fatigue
- Ejector with infinitely variable length adjustment makes it easier for different tips
- Low pipette weight
- Soft springs system provides comfort at work
- Fully autoclavable
- UV resistant
- Calibration key included
- Factory-calibrated to strict ISO 8655 standards

Single channel pipettors

Cat. No.	Volume range (µl)
12740	0,5-10
12750	2-20
12770	10-100
12780	20-200
12790	100-1000

8-channel pipettors

Cat. No.	Volume range (µl)
12800	1-10
12820	5-50
12460	20-200
12840	50-300

