

Distributed by



Made of natural virgin  
polypropylene



Universal tips



Graduated tips



RNase/DNase free



No mould release  
agent used




Bi-component tips: the upper  
part is a flexible elastomer for  
optimum sealing and precision

Tips ideal for multi-channel



Pointes Multi Sorenson



Recycling  
process



Reduction of  
muscle fatigue

Each tip is housed in an individually sealed cell,  
eliminating any risk of tip-to-tip contamination.  
Available with or without filter

## CellRack

- Innovative : Bi-component tips : Flexible elastomer and polypropylene
- Hold tips steadily in fully enclosed cells, minimize static and contamination
- Excellent platforms for single and multichannel pipetting
- Recution of muscle fatigue : less pressure to create the optimal seal and less force required for tip ejection



### Sterile filter tips

Volume (µl)	Cat.No.	Units/box
1 - 20	010320K	5 x 960
10 - 100	010510K	5 x 960
20 - 330	010350K	5 x 960
100 - 1250	010380k	5 X 768

### Sterile non-filter tips

Volume (µl)	Cat.No.	Units/box
1-20	010310K	5 x 960
20 - 330	010340K	5 x 960
100 - 1250	010370k	5 X 768

### Non-Sterile non-filter tips

Volume (µl)	Cat.No.	Units/box
1-20	010300K	5 x 960
20 - 330	010330K	5 x 960
100 - 1250	010360K	5 X 768