




**BETTER BASICS**  
LABORBEDARF

Distributed by

**BIOMEDICA**  
www.bmgrp.eu



1<sup>ST</sup> QUARTER 2025



# Product catalogue

SmartRack®

Unlock Your Lab Potential

# Contents

Your Advantage in Everyday Laboratory Work.....	3
SmartRack®   Mobile Organizationsystem.....	6
The Configuration.....	6
Pre-configurations.....	7
Frame size.....	17
Additional accessories.....	19
SmartCase.....	22
Customizing.....	23
Item overview.....	25
Technical information.....	26

# Your Advantage in Everyday Laboratory Work

## Better Basics Laborbedarf – Your Advantage in Everyday Laboratory Work

Better Basics Laborbedarf optimises laboratories and their workspaces with innovative, award-winning organisational systems. We support our customers in creating the best possible working environment for their laboratory staff, helping to unlock the full potential of their laboratories. This approach ensures consistently high quality, efficiency, and process stability in the lab while providing a safe, ergonomic workplace with satisfied employees.

### Enhancing Employee Satisfaction

Our SmartRack® product familie offer a comprehensive system to optimise laboratory workflows. With over 100 modular solutions, they not only ensure order but also noticeably improve working conditions. Enhanced ergonomics, reduced movement paths, and clearly structured workstations increase both efficiency and employee satisfaction. Safe and organised environments reduce stress and promote job satisfaction. Furthermore, they demonstrably contribute to improved workplace safety and occupational health through thoughtful placement and ease of use.


### A System for Sustainable Process Improvement

Our product families are designed to standardise workflows and save time. The modular concept of the SmartRack® allows for flexible adaptations to any laboratory environment – from safety workbenches to large open-plan labs. These systems not only optimise space utilisation but also increase efficiency by an average of 30% and save up to 80% space. With a wide range of modules and accessories, they offer customised solutions for every need.

### Efficient, Safe, and Satisfying

For decision-makers and managers, the SmartRack® ensures process stability and economical resource use. At the same time, these systems provide innovation managers and organisational developers with the ideal basis for forward-looking workplace designs. For laboratory staff, this means fewer errors thanks to well-structured work areas, a more relaxed working environment, and often the ability to finish earlier due to significantly improved workflows.

**Discover how Better Basics Laborbedarf can transform your everyday laboratory work – for more efficient processes, satisfied employees, and excellent results.**



### ▶ SmartRack®, the mobile solution for: laboratory workbench, fume hood, safety cabinet



### + Customised solutions

### ▶ Modules

- Currently, more than 100 modules
- All modules in the Better Basics catalogue can be used universally in all SmartRack® solutions
- Customized modules can be developed based on the specific customer needs, all customized modules can be used in all SmartRack® solutions



### + Customised solutions

▶ **SmartRack®**, the mobile solution for: laboratory workbench, fume hood, safety cabinet



✚ Customised solutions

## The SmartRack® – Efficiency and Safety Conveniently Combined.

The SmartRack® is the universal solution for laboratory organisation. With its modular and flexible structure, it optimises workspace utilisation and achieves up to 80% space savings. Whether in a laboratory fume hood, on the workbench, or in a biosafety cabinet – the SmartRack® adapts to your needs.

It supports efficient workflows, minimises sources of error, and enhances safety. Its pre-configured variants and customisable elements make the SmartRack® the ideal foundation for any laboratory environment.



Photo: A SmartRack® classic inside a WALDNER SCALA laboratory fume hood

# SmartRack® | The Configuration

## 1. Selection of pre-configurations

Basic Equipment, Biology, Chemistry, Pipetting Set, Consumables, etc.



OR

## 2. Customised configuration

### 2.1 Frame size and frame type selection:

SmartRack® classic (mini) or SmartRack® primus (mini)



### 2.2. Selection of individual modules:

holders and storage box and storage box inserts



### 2.3 Selection of additional SmartRack® accessories (optional):

if required further slot rail, drip tray or protective cover



## Configuration of a personalized SmartRack®, SmartRack® classic (left) plus modules (right)



Image above: SmartRack® classic



Images above: Selection of modules

**Each SmartRack® classic** (article number 01704) or **SmartRack® classic mini** (article number 01707) **is already delivered ex works with the following accessories**

- 1 x Pair of feet
- 1 x Drip tray
- 3 x Slot rails each with handles
- 1 x Frame

## SmartRack® mini Pipetting Set



The SmartRack® mini Pipetting Set has been specially designed for working with liquids in small quantities and includes the following components (item number of the individual components in brackets for better understanding):

9 components:

- 1 x SmartRack® classic mini (Item no. 01707)
- 1 x Holder for bin liner (Item no. 01202)
- 1 x Holder for pipettes - for 3 pipettes with hook (Item no. 01145)
- 2 x Holder for pipette tip boxes (Item no. 01147)
- 2 x Holder for laboratory screw top bottles, 100 - 1000 ml (GL45) (Item no. 01122)
- 2 x Holder for well plates (Item no. 01125)

Detailed views of the individual components included:



Item no. 01707



Item no. 01202



Item no. 01145



Item no. 01147



Item no. 01122



Item no. 01125

### Note on the article numbers provided:

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

Trash bag, pipettes, pipette tip boxes, laboratory screw top bottle and well plates not included in scope of delivery

Item	Type	PU	Item no.
SmartRack® mini Pipetting Set, bundle price:		1	01512

## SmartRack® Basic Equipment



The SmartRack® Basic Equipment is a basic set for organising basic accessories in everyday laboratory use and includes the following components (item number of the individual items in brackets for better understanding):

6 components:

- 1 x SmartRack® classic (Item no. 01704)
- 1 x Holder for bin liner (Item no. 01202)
- 1 x Storage box (Item no. 01142)
- 1 x Holder for kitchen roll (Item no. 01162)
- 1 x Holder for glove / tissue boxes (Item no. 01143)
- 1 x Holder for wash bottles, for 500 ml (up to Ø 75 mm) (Item no. 01146)

Detailed views of the individual components included:



Item no. 01704



Item no. 01202



Item no. 01142



Item no. 01162



Item no. 01143



Item no. 01146

### Note on the article numbers provided:

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

Small accessories, trash bag, tissue box, wash bottle and kitchen roll not included in scope of delivery

Item	Type	PU	Item no.
SmartRack® Basic Equipment, bundle price:		1	01506



## SmartRack® Starter Set Biology



The SmartRack® Starter Set Biology has been specially designed for work in the field of biology and includes the following individual components (item number of the individual component in brackets for better understanding):

8 components:

- 1 x SmartRack® classic (Item no. 01704)
- 2 x Holder for laboratory screw top bottles, 100 - 1000 ml (GL 45) (Item no. 01122)
- 2 x Storage box (Item no. 01142)
- 1 x Holder for pipettes - for 3 pipettes with hook (Item no. 01145)
- 1 x Holder for pipette tip boxes (Item no. 01147)
- 1 x Holder for microliter syringes (Item no. 01155)

Detailed views of the individual components included:



Item no. 01704



Item no. 01122



Item no. 01142



Item no. 01145



Item no. 01147



Item no. 01155

### Note on the article numbers provided:

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

Laboratory screw top bottle, pipettes, pipette tip box and microliter syringes not included in scope of delivery

Item	Type	PU	Item no.
SmartRack® Starter Set Biology, bundle price:		1	01507

# SmartRack® | Pre-configurations

## SmartRack® Starter Set Chemistry



The SmartRack® Starter Set Chemistry has been specially designed for work in the field of chemistry and includes the following components (item number of the individual components in brackets for better understanding):

8 components:

- 1 x SmartRack® classic (Item no. 01704)
- 2 x Holder for glassware with ST14, plus closure (Item no. 01112)
- 2 x Holder for glassware with ST29, plus closure (Item no. 01115)
- 2 x Holder for chromatography column, Ø 25 - 35 mm (Item no. 01107)
- 1 x Storage box (Item no. 01142)

Detailed views of the individual components included:



Item no. 01704



Item no. 01112



Item no. 01115



Item no. 01107



Item no. 01142

### **Note on the article numbers provided:**

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

Glassware and chromatography column not included in scope of delivery

Item	Type	PU	Item no.
SmartRack® Starter Set Chemistry, bundle price:		1	01508

## SmartRack® for Consumables



The SmartRack® for Consumables has been specially designed for the provision of consumables at the laboratory workplace and comprises the following components (item number of the individual components in brackets for better understanding):

12 components:

- 1 x SmartRack® classic (Item no. 01704)
- 1 x Dispenser for Parafilm® rolls (Item no. 01171)
- 1 x Dispenser for adhesive tapes (Item no. 01172)
- 3 x Holder for glove / tissue boxes (Item no. 01143)
- 1 x Holder for bin liner (Item no. 01202)
- 3 x Holder for wash bottles, for 500 ml (up to Ø 75 mm) (Item no. 01146)
- 1 x Holder for kitchen roll (Item no. 01162)
- 1 x Slot rail with handles (Item no. 01601)



# SmartRack® | Pre-configurations

Detailed views of the individual components included:



Item no. 01704



Item no. 01171



Item no. 01172



Item no. 01143



Item no. 01202



Item no. 01146



Item no. 01162



Item no. 01601

**Note on the article numbers provided:**

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

*Parafilm® rolls, adhesive tapes, glove / tissue boxes, bin liner, wash bottles and kitchen roll not included in scope of delivery*

Item	Type	PU	Item no.
SmartRack® for Consumables, bundle price:		1	01509

## SmartRack® for Life Science



The SmartRack® for Life Science has been specially designed for work in the fields of biology, biomedicine, biotechnology as well as special fields of biology and for handling liquids in small quantities and preparing and storing samples for later analysis. It includes the following components (item number of each component in brackets for better understanding):

16 components:

- 1 x SmartRack® classic (Item no. 01704)
- 1 x Holder for bin liner (Item no. 01202)
- 1 x Storage box (Item no. 01142)
- 1 x Storage box insert - for small accessories (Item No. 01148)
- 1 x Storage box insert - for microreaction vessels, 20 x 2 ml (Item no. 01130)
- 1 x Storage box insert - for centrifuge tubes, 6 x 15 ml, 4 x 50 ml (Item no. 01102)
- 2 x Holder for laboratory screw top bottles, 100 - 1000 ml (GL45) (Item no. 01122)
- 2 x Holder for well plates (Item no. 01125)
- 1 x Holder for wash bottles, for 500 ml (up to Ø 75 mm) (Item no. 01146)
- 1 x Holder for pipettes - for 3 pipettes with hook (Item no. 01145)
- 2 x Holder for pipette tip boxes (Item no. 01147)
- 1 x Dispenser for Parafilm® rolls (Item no. 01171)
- 1 x Slot rail with handles (Item no. 01601)



# SmartRack® | Pre-configurations

Detailed views of the individual components included:



Item no. 01704



Item no. 01202



Item no. 01142



Item no. 01148



Item no. 01122



Item no. 01125



Item no. 01146



Item no. 01145



Item no. 01147



Item no. 01171



Item no. 01130



Item no. 01102



Item no. 01601



**Note on the article numbers provided:**

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

Bin liner, small accessories, laboratory screw top bottles, well plates, wash bottle, pipettes, pipette tip boxes, parafilm® roll, microreaction vessels and centrifuge tubes not included in scope of delivery

Item	Type	PU	Item no.
SmartRack® for Life Science, bundle price:		1	01510

## SmartRack® for Chromatography



The SmartRack® for Chromatography has been specially designed for substance separation tasks by means of column chromatography and comprises the following components (item number of the individual components in brackets for better understanding):

17 components:

- 1 x SmartRack® classic (Item no. 01704)
- 1 x Storage box (Item no. 01142)
- 1 x Storage box insert – for small accessories (Item no. 01148)
- 1 x Storage box insert - for sample containers and test tubes, 12x Ø 10 - 20 mm (Item no. 01160)
- 1 x Storage box insert - for sample containers and test tubes, 20x Ø 3 - 10 mm (Item no. 01161)
- 6 x Holder for chromatography column, Ø 25 - 35 mm (Item no. 01107)
- 1 x Holder for wash bottles, for 500 ml (up to Ø 75 mm) (Item no. 01146)
- 1 x Holder for microliter syringes (Item no. 01155)
- 2 x Holder for laboratory screw top bottles, 100 - 1000 ml (GL45) (Item no. 01122)
- 1 x Holder for glove / tissue boxes (Item no. 01143)
- 1 x Slot rail with handles (Item no. 01601)



# SmartRack® | Pre-configurations

Detailed views of the individual components included:



Item no. 01704



Item no. 01142



Item no. 01107



Item no. 01146



Item no. 01155



Item no. 01122



Item no. 01148



Item no. 01160



Item no. 01161



Item no. 01143



Item no. 01601

**Note on the article numbers provided:**

The modules listed here can be found in the Better Basics partial catalog "Modules." You can download this catalog at [www.download-en.better-basics-laborbedarf.de](http://www.download-en.better-basics-laborbedarf.de).

Chromatography columns, wash bottle, microliter syringes, laboratory screw top bottles, small accessories, sample containers and test tubes, glove / tissue box not included in scope of delivery

Item	Type	PU	Item no.
SmartRack® for Chromatography, bundle price:		1	01511



## SmartRack® classic



The SmartRack® classic was developed for use at many different laboratory workstations and especially for work in fume hoods. The SmartRack® classic mini is specially designed for work in safety cabinets. It consists of an aluminium frame with a black epoxy-polyester-based powder coating. Every third slot rail holder on both sides of the SmartRack® classic frame has a recess at the front to make it easier to align the slotted rails horizontally.

SmartRack® classic or SmartRack® classic mini includes the following components:

- 1 x frame incl. one pair of rear shortenable feet with tray holder
- 3 x slot rails each with handles
- 1 x drip tray

Size comparison SmartRack® classic (left) and SmartRack® classic mini (right):



Item	Type	PU	Item no.
SmartRack® classic	635 x 615 x 265 mm (WxHxD)	1	01704
SmartRack® classic mini	420 x 380 x 265 mm (WxHxD)	1	01707

The specified width of the SmartRack® and SmartRack® mini includes the frame, the slot rails with handles, as well as the feet.

If you need a SmartRack® with special size requirements or custom brackets, please contact us at [contact@better-basics.de](mailto:contact@better-basics.de).

# SmartRack® | Frame size and frame type

## SmartRack® primus



The SmartRack® primus is the affordable entry-level version into the smart world of laboratory equipment from Better Basics. It was developed for use at many different laboratory workstations and in particular for work at the fume hood. The SmartRack® primus mini is specially designed for working on the safety cabinet. It consists of an uncoated, polished aluminium frame. Every third slot rail holder on both sides of the SmartRack® primus frame has a recess at the front to make it easier to align the slotted rails horizontally.

The SmartRack® primus or SmartRack® primus mini comprises the following components:

- 1 x frame incl. one pair of aluminium feet
- 2 x slotted rails incl. two handles on each side

Size comparison SmartRack® primus (left) and SmartRack® primus mini (right):



Item	Type	PU	Item no.
SmartRack® primus	635 x 615 x 265 mm (WxHxD)	1	01711
SmartRack® primus mini	420 x 380 x 265 mm (WxHxD)	1	01712

The specified width of the SmartRack® and SmartRack® mini includes the frame, the slot rails with handles, as well as the feet.

If you need a SmartRack® with special size requirements or custom brackets, please contact us at [contact@better-basics.de](mailto:contact@better-basics.de).

## Drip tray



The drip tray prevents the spread of liquids. In case of contamination, the drip tray can be easily removed from the SmartRack® and cleaned.

Item	Type	PU	Item no.
Drip tray	for SmartRack®	1	01603
Drip tray mini	for SmartRack® mini	1	01604


## Slot rail with handles



The slot rail can be used in many different positions in the SmartRack®. A larger number of slot rails offers more options for placing modules through an optimal vertical arrangement, providing more scope for customizing the structure of the SmartRack® and adapting the module arrangement to the individual needs of each user.

Item	Type	PU	Item no.
Slot rail with handles	for SmartRack®	1	01601
Slot rail mini with handles	for SmartRack® mini	1	01602

## Back-to-back connectors



The back-to-back connectors for the SmartRack® enable two SmartRack® to be securely connected and thus used as a isolated solution. With the rear foot sections removed, two SmartRack® can be securely connected by simply attaching the connectors to the horizontal or diagonal braces of the frame rear panel. A set consists of three individual back-to-back connectors.

Item	Type	PU	Item no.
Back-to-back connectors for SmartRack®		1	01614

## Protective cover for SmartRack®



Version with zippers

The protective hood made of fiber-reinforced polyethylene (mesh film), which was specially designed for the SmartRack®, reliably protects the SmartRack® against dirt and dust particles as well as providing splash protection against aggressive liquids. The dust and particle protection hood is available in two versions that do not differ in size. As a basic model and as a model with two zippers on the edges of the front. For regular removal of the dust and particle protection hood and for use in unfavorable height conditions, we recommend choosing the model with zippers.

*SmartRack® not included in scope of delivery*



Version without zippers



SmartRack® mini cover with zippers

Item	Type	PU	Item no.
Protective cover for SmartRack®	with zipper	1	01609
Protective cover for SmartRack®	without zipper	1	01608
Protective cover for SmartRack® mini	with zipper	1	01613
Protective cover for SmartRack® mini	without zipper	1	01612

# SmartCase

## SmartCase



SmartRack® not included, other equipment not included

The SmartCase is designed to transport and store the SmartRack®, a modular organization system for laboratory workstations. Its aluminium case safely protects a SmartRack® inside and its contents from light and other environmental influences. The SmartCase consists of a lower part with four castors into which a SmartRack® can be placed. The upper part of the SmartCase has a cushion on the inside, which serves to fix the SmartRack® in the SmartCase as soon as the SmartCase is closed by the tension locks. This prevents the SmartRack® from slipping or tipping over. The SmartCase has a self-acting valve system which can be used to change the atmosphere inside the SmartCase by means of an external gas supply. This allows, for example, the temporary generation of a low-oxygen atmosphere in the SmartCase, which creates an additional safety level when transporting oxygen-sensitive, highly flammable substances or substances that spontaneously ignite in air. For further information on the durability of the atmosphere exchange in the SmartCase, please contact us at [contact@better-basics.de](mailto:contact@better-basics.de).

The SmartCase includes the following components:

- A bottom section with four casters, four tension locks and an automatic gas outlet valve system.
- A top part with a handle, four put-on aids and an automatic gas inlet valve system
- A connector for connecting gas hoses of various diameters to the SmartCase

Item	Type	PU	Item no.
SmartCase		1	02101

## Customization of SmartRack® Products

Better Basics Laborbedarf offers an unparalleled approach to laboratory workstation organisation with its customisable and modular solutions. Below is an overview of how our customisation services span the SmartRack® product line:

### 1. Customization of the SmartRack®

#### a) For the SmartRack® and the SmartRack® mini

Our SmartRack® system is fully customisable to suit the specific needs of your workspace. Both the standard and mini versions can be adapted to the dimensions of your equipment and workflow. From compact configurations to larger frameworks, we enhance your laboratory's organisation and efficiency.

#### b) Frame Size Adjustments

Adjustable widths and heights allow the SmartRack® to integrate seamlessly into fume hoods, safety workbenches, or experimental setups. Each tailored frame is paired with appropriately sized slot rails, ensuring modular compatibility and reliability.

### 2. Custom Modules

#### a) Tailored for Any Laboratory Equipment

Using advanced 3D printing technology, we create modules and holders in any shape or size to fit your laboratory equipment. From pipette holders to mounts for specialised glassware, our bespoke designs cater to your unique requirements. Our rapid prototyping process ensures swift delivery of tailored solutions.

### Why Choose Better Basics for Customisation?

- 1. Precision Craftsmanship:** Each component is tailored to your exact specifications, guaranteeing optimal performance and seamless integration.
- 2. Fast Delivery:** Our in-house design and manufacturing teams ensure swift turnaround from concept to delivery.
- 3. Transparent Pricing:** Costs for customisation are agreed upon in advance, providing clarity and control.
- 4. Scalability:** Whether expanding an existing laboratory or building a new one, our products adapt to meet your growing needs.

### Transform Your Laboratory Workspace

With Better Basics Laborbedarf, your laboratory can benefit from customised SmartRack® products. Visit [customising.better-basics-laborbedarf.de](https://www.better-basics-laborbedarf.de) or contact our team today to unlock the full potential of your laboratory space.

This concise yet comprehensive approach ensures your laboratory benefits from enhanced efficiency, organisation, and safety—all tailored to your specific needs.

# Customizing

## Examples for realised customisations based on customer requests

### Special frame sizes

#### SmartRack® frame - customisable in width as well as height



#### SmartRack® frame - extra wide

For a special requirement we manufactured an extra wide frame of the SmartRack®. Due to the width, this SmartRack® makes optimal use of the working space in our customer's laboratory.

### Special modules for laboratory accessories and laboratory glassware

#### Holder for special glassware



#### Double holder for glassware with ST10

For a project, a customer needed a double holder for glassware ST10. The module allows glassware with standard ground joint ST10 to be easily hung next to each other. The specially shaped recess prevents the glassware from falling out.

Of course, adaptations of this kind are also conceivable for other glassware.

### Special modules for laboratory equipment

#### Holder for dandyVice wall mount



In cooperation with dandy innovations, a module for mounting a dandyVice wall bracket was developed. The combination of dandyVice wall mount and holder module from Better Basics Laborbedarf makes it easy to use the dandyVice in the SmartRack®.



## SmartRack® | Pre-configurations

Item	Type	Item no.	PU	Page
SmartRack® Basic Equipment		01506	1	8
SmartRack® Starter Set Biology		01507	1	9
SmartRack® Starter Set Chemistry		01508	1	10
SmartRack® for Chromatography		01511	1	15
SmartRack® for Consumables		01509	1	11
SmartRack® for Life Science		01510	1	13
SmartRack® mini Pipetting Set		01512	1	7

## SmartRack® | Frame

Item	Type	Item no.	PU	Page
SmartRack® classic	635 x 615 x 265 mm (WxHxD)	01704	1	17
SmartRack® classic mini	420 x 380 x 265 mm (WxHxD)	01707	1	17
SmartRack® primus	635 x 615 x 265 mm (WxHxD)	01711	1	18
SmartRack® primus mini	420 x 380 x 265 mm (WxHxD)	01712	1	18

## SmartRack® | Additional accessories

Item	Type	Item no.	PU	Page
Back-to-back connectors for SmartRack®		01614	1	20
Drip tray	for SmartRack®	01603	1	19
Drip tray mini	for SmartRack® mini	01604	1	19
Protective cover for SmartRack® mini	with zipper	01613	1	21
Protective cover for SmartRack® mini	without zipper	01612	1	21
Protective cover for SmartRack®	with zipper	01609	1	21
Protective cover for SmartRack®	without zipper	01608	1	21
Slot rail mini with handles	for SmartRack® mini	01602	1	19
Slot rail with handles	for SmartRack®	01601	1	19

## SmartCase

Item	Type	Item no.	PU	Seite
SmartCase		02101	1	22

Net prices plus statutory VAT. | Valid from 01.01.2025

The general terms and conditions of Better Basics Laborbedarf GmbH apply.

# Technical information

## Technical information:

- All products are washable with water and mild detergent. In addition all components can be wiped with alcohol-based disinfectants.
- The plastics and metals used are 100% recyclable.
- Our products are autoclavable (maximum temperature 121 °C).
- The plastic elements are resistant to heat deformation under normal laboratory conditions.
- To gently clean the anodised aluminium rail, we recommend using a pH-neutral hand-washing detergent. Ensure that the hand wash detergent is odourless and contains no perfumes, greases or additives. Clean the slot rail by hand and rinse thoroughly with clean water.

## Note on joint tolerances for ST14:

- Maximum outside diameter at the narrowest point of the cone joint: 17.2 mm (larger diameters are too wide for the holder and cannot be held)
- Minimum outer diameter at the widest point of the cone joint: 18.0 mm (smaller diameters are too narrow for the holder and would fall out)

## Note on joint tolerances for ST29:

- Maximum outside diameter at the narrowest point of the cone joint: 31.4 mm (larger diameters are too wide for the holder and cannot be held)
- Minimum outer diameter at the widest point of the cone joint: 32.8 mm (smaller diameters are too narrow for the holder and would fall out)

**IMPORTANT TO NOTE:** The SmartRack® and SmartRack® mini as well as the SmartCase and SmartStand are not suitable for the permanent storage of hazardous substances.

We reserve the right to make technical changes, material or model changes, errors and misprints excepted. If you have any further questions, please do not hesitate to contact us at [contact@better-basics.de](mailto:contact@better-basics.de)

## For more detailed information, please visit our website:

Functionality of SLS 3D printing:	<a href="http://www.sls-3d-printing.better-basics-laborbedarf.de">www.sls-3d-printing.better-basics-laborbedarf.de</a>
Resistance to UV light:	<a href="http://www.uv-resistance.better-basics-laborbedarf.de">www.uv-resistance.better-basics-laborbedarf.de</a>
Resistance to aggressive media:	<a href="http://www.material-resistance.better-basics-laborbedarf.de">www.material-resistance.better-basics-laborbedarf.de</a>
Sustainability and design-efficient manufacturing:	<a href="http://www.sustainability.better-basics-laborbedarf.de">www.sustainability.better-basics-laborbedarf.de</a>



Image: Building of the Better Basics Laborbedarf GmbH in Dresden

## YOU ARE INTERESTED IN OUR PRODUCTS OR YOU HAVE FURTHER QUESTION?

Get in touch with us:

[contact@better-basics.de](mailto:contact@better-basics.de)

You can also find more information and news about our product innovations at

[www.Better-Basics-Laborbedarf.de](http://www.Better-Basics-Laborbedarf.de)

### **Better Basics Laborbedarf GmbH**

Loebtauer Str. 69

01159 Dresden

Germany

T: +49 (0) 178 9617 577 (Direct customers)

T: +49 (0) 176 6233 8026 (Sales partners)

E: [contact@better-basics.de](mailto:contact@better-basics.de)

W: [www.Better-Basics-Laborbedarf.de](http://www.Better-Basics-Laborbedarf.de)

### **Our request to you**

Our products such as our SmartRack® are only as good as their users. With our roots in the Leibniz Institute for Polymer Research Dresden, our company Better Basics Laborbedarf stands for the tradition of "Made in Germany" quality. We feel in multiple aspects committed to this quality label for German companies, which stands for high innovative strength. Even beyond our product catalogue, the same applies for us at Better Basics: Just contact us if you have ideas for extensions or improvements or if you notice something in your daily laboratory routine that we should optimize right away.

We want to create the best possible working environment for you in the laboratory and we are grateful for any comments and inspiration. We can jointly shape the future of laboratory work by cooperating with scientists and researchers. Just like you, we're always looking for a way to make our world a little better every day.

### **Impress:**

Better Basics Laborbedarf GmbH  
Loebtauer Str. 69, G-01159 Dresden

E-Mail: [info@better-basics.de](mailto:info@better-basics.de)

CEO: Mario Schneider

Commercial register: HRB 39524; Register court: Amtsgericht Dresden

VAT registration number: DE327508812

Distributed by

