



**DOMINIQUE
DUTSCHER
DISTRIBUTION**

Distributed by

BIOMEDICA
www.bmgrp.eu



15 and 50 ml centrifuge tubes

- Polypropylene tubes and polyethylene caps
- Standard and reverse graduation
- Graduation on the bottom of the tube
- Indication of maximum freezing level and maximum centrifuging speed
- Large marking area

Volume (ml)	15		50		
Ext diameter (mm)	15.4		27.8		
Height (mm)	119		115		
Autoclavable	Yes				
Max centrifugal speed (g)	15000		12000		
Skirt	No	No	No	No	Yes
Colour	Transparent				
Operating temperature	- 86 °C to + 121 °C				
Specifications	Free of DNAs, RNAs and endotoxins; perfectly watertight at - 95 kPa				
Units/box	20 bags of 25 tubes	10 racks of 50 tubes	20 bags of 25 tubes	20 racks of 25 tubes	20 bags of 25 tubes
Part no. turquoise cap	923015	923016	923012	923013	923014
€/box	NC -	NC -	NC -	NC -	NC -



Serological pipettes

Serological pipettes

- In high transparency polystyrene
- Easy-to-read bi-directional scale
- Sterile (SAL 10⁻⁶)
- Colour-coded on pack and pipette
- With synthetic filter
- Anti-drip system on 25 and 50 ml models
- Available in bulk or individually wrapped
- Part number, expiry date and batch number indicated on individual packaging
- 1, 2, 5 and 10 ml models in one piece



Vol. (ml)	Graduation (ml)	Colour code	Individual paper/plastic packaging				In bulk			
			Part no.	Units/bag	Units/box	€/ box	Part no.	Units/bag	Units/box	€/ box
1	0.01	Yellow	923000	100	1000	NC -	923001	25	1000	NC -
2	0.01	Green	923002	100	800	NC -	923003	25	800	NC -
5	0.1	Blue	923004	50	300	NC -	923005	25	500	NC -
10	0.1	Orange	923010	50	200	NC -	923011	25	350	NC -
25	0.2	Red	923006	50	200	NC -	923007	25	350	NC -
50	1	Purple	923008	25	100	NC -	923009	25	100	NC -
100	1	Purple	923033	-	50	NC -	-	-	-	-



Short serological pipettes

- Length 232 mm
- Polystyrene
- With cotton
- Colour code by volume
- Bidirectional dimming
- Easy to handle under a fume hood
- Precision identical to that of a standard pipette



Vol. (ml)	Graduation (ml)	Colour code	Individual paper/plastic packaging			In a bag		
			Part no.	Units/box	€/ box	Part no.	Units/box	€/ box
5	0.1	Blue	923017	200	NC -	923020	200	NC -
10	0.2	Orange	923018	200	NC -	923021	200	NC -
25	0.5	Red	923019	200	NC -	923022	100	NC -
50	1	Purple	-	-	-	923023	50	NC -

